



UAB "LT Biotech"

J.K. 302303586, PVM kodas LT 100004741118,
Rugiu 21-24, LT-08419, Vilnius Reg. Nr. 127918, V.I.
Registru centras Vilniaus filialas
Tel/fax +370 5 216 02 27

Hanks'BalancedSalts,w/oCa&Mg,w/oPhenolRedLTBSS02

General Information

Hank'sBalancedSaltSolution(HBSS)isusedforavarietyofcellcultureapplications,suchas rinsingcells,transportingcellsortissuesamples, dilutingcellsforcounting,andpreparing reagents.TheessentialfunctionofabalancedsaltsolutionistomaintainpHandosmoticbalanceas well as providing cells with water and essential inorganic ions.

Product Specifications

Appearance	Clearliquid
Storage and shelf life	Storeat+2to+25°Cuponarrival. Onceopenedstoreat+2to+25°Candusewithin24month.
Shipping conditions	Ambient

Formulation

Component	concentration mg/l
D-Glucose	1000
KCl	400
KH ₂ PO ₄	60
NaCl	8000
NaHCO ₃	350
Na ₂ HPO ₄	48

Precautions and Disclaimer

Thisproductisforresearchuseonly.